

## Congratulations on the purchase of your Planmeca Compact i Classic Treatment Unit.

Henry Schein will be your point of contact for any parts or consumables you may require, and the list provided below can be used to help identify specific products for routine maintenance, cleaning and service parts for your new Planmeca chair. All your parts and consumable items can either be ordered via the Henry Schein Customer Services Team or through the **Henry Schein Webshop - [www.henryschein.co.nz](http://www.henryschein.co.nz)**.

### Planmeca Compact i Classic Treatment Unit Accessories

#### Suction Filter (Strainer)

PB-00221013 - Pack of 1



#### Luzzani Syringe Tip

PB-00871348

#### Luzzani Metallic Cone

PB-00871441



#### Luzzani Minibrite Syringe Tip

PB-10037530



#### 16.5mm High Volume Suction Tubing

PB-02200192 - 2.0m Grey (Tubing Only)

#### 10.0mm Saliva Ejector Suction Tubing

PB-02200191 - 2.0m Grey (Tubing Only)



#### High Volume Suction Adaptor

PB-00251002 - 16mm > 11mm Reduction



#### Light Handle for Solanna Light

PB-FE004325 - Removable, Sterilisable Each

